


Shuttle-Flex



Features

- + Tension bands provide a **dampening** function at heel strike and smooth heel rocker. During standing, they provide a **centering**, stabilizing force



Segment attaches to the Shuttle ankle joint stirrup

Standard

Available **separately or as combi-product** together with the Shuttle ankle joint

Precise assembly calibrated to 1/100mm for smooth operation free of play


Serial number for traceability

Clearly illustrated instructions for use

Fig. left













We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product		AFO/KAFO	HKAFO/RGO		Article No.
Shuttle-Flex-22		up to 15 kg	up to 15 kg		SH-Flex-22
Shuttle-Flex-25		up to 25 kg	up to 25 kg		SH-Flex-25
Shuttle-Flex-28		up to 35 kg	up to 35 kg		SH-Flex-28
Shuttle-Flex-30		up to 45 kg	up to 45 kg		SH-Flex-30
Shuttle-Flex-32		up to 55 kg	up to 55 kg		SH-Flex-32
Shuttle-Flex-34		(up to 70 kg) not yet available			SH-Flex-34

									
Adjustable Alignment	Modular	Force Options	Range of Motion	Swing Phase	Quick Adjustment	Interchangeable Stops	No Grinding Necessary	Noise Dampening	Water Resistant