

Shuttle-Uno-back




Features

Posterior power unit uses medium to strong forces (plantarflexion resistance) for **reducing** correctable equinus. Weaker forces are used to dampen heel strike for **a smooth heel rocker**

Tension bands provide **a smooth rollover** during stance phase

Segment attaches to the Shuttle ankle joint stirrup



Standard


Available **separately or as combi-product** together with the Shuttle ankle joint

Precise assembly calibrated to 1/100mm for smooth operation free of play


Serial number for traceability

Clearly illustrated instructions for use

Fig. left











We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

| Article No. Combi-product |  | AFO/KAFO |  | Article No. |
|---------------------------|---|---------------------------------|---|----------------|
| Shuttle-Uno-back-22 | | up to 15 kg | | SH-Uno-back-22 |
| Shuttle-Uno-back-25 | | up to 25 kg | | SH-Uno-back-25 |
| Shuttle-Uno-back-28 | | up to 35 kg | | SH-Uno-back-28 |
| Shuttle-Uno-back-30 | | up to 45 kg | | SH-Uno-back-30 |
| Shuttle-Uno-back-32 | | up to 55 kg | | SH-Uno-back-32 |
| Shuttle-Uno-back-34 | | (up to 70 kg) not yet available | | SH-Uno-back-34 |

| | | | | | | | | | |
|---|---|---|---|--|------------------|---|---|---|-----------------|
|  |  |  |  |  | Quick Adjustment |  |  |  | Water Resistant |
| Adjustable Alignment | Modular | Force Options | Range of Motion | Swing Phase | | Interchangeable Stops | No Grinding Necessary | Noise Dampening | |