

# Liber-Extensor-N



## Feature

+ **Force unit** assists knee extension and helps when sitting down and standing up


+ **2-piece adjustable body** to vary flexion stop between **95°** and **125°**

+ **7 positions** to shift lever arm and increase or decrease torque

**Cushioned** flexion stop

Wide range of **force options** available

**2 different styles** of base plate for versatility: retrofit or incorporated



## Standard


**Fully assembled** complete with everything needed for a functional add-on to the Liber knee joint

Available separately or as **combi-product** together with the Liber knee joint


**Serial number** for traceability

**Clearly illustrated** instructions for use



Fig. left



We recommend using our configurator to make your selection.  
The configurator considers, among other factors, weight to height ratios.



#### Weight Classification (body weight)

Article No. Combi-product		KAFO	HKAFO		Article No.
Liber-Extensor-N-22		up to 15 kg	up to 15 kg		Extensor-N-22
Liber-Extensor-N-25-Plus		up to 30 kg	up to 25 kg		Extensor-N-25-Plus
Liber-Extensor-N-28		up to 40 kg	up to 35 kg		Extensor-N-28
Liber-Extensor-N-30		up to 50 kg	up to 45 kg		Extensor-N-30
Liber-Extensor-N-30-Plus		up to 70 kg	up to 60 kg		Extensor-N-30-Plus
Liber-Extensor-N-34		up to 85 kg	up to 85 kg		Extensor-N-34

								
Locked	Extension	Flexion	Posterior Offset	Modular	Force Options	No Grinding Necessary	Noise Dampening	Water Resistant