

Liber-Extensor-S



Features

+ **Force unit** assists knee extension and helps when sitting down and standing up


+ **3 positions** to shift lever arm and increase or decrease torque

+ **Flexion up to 115°** for comfortable sitting

Cushioned flexion stop

Wide range of **force options** available

2 different styles of base plate for versatility: retrofit or incorporated



Standard


Fully assembled complete with everything needed for a functional add-on to the Liber knee joint

Available separately or as **combi-product** together with the Liber knee joint


Serial number for traceability

Clearly illustrated instructions for use



Fig. left



We recommend using our configurator to make your selection.
The configurator considers, among other factors, weight to height ratios.



Weight Classification (body weight)

Article No. Combi-product		KAFO	HKAFO		Article No.
Liber-Extensor-S-22		up to 15 kg	up to 15 kg		Extensor-S-22
Liber-Extensor-S-25-Plus		up to 30 kg	up to 25 kg		Extensor-S-25-Plus
Liber-Extensor-S-28		up to 40 kg	up to 35 kg		Extensor-S-28
Liber-Extensor-S-30		up to 50 kg	up to 45 kg		Extensor-S-30
Liber-Extensor-S-30-Plus		up to 70 kg	up to 60 kg		Extensor-S-30-Plus
Liber-Extensor-S-34		up to 85 kg	up to 85 kg		Extensor-S-34

								
Locked	Extension	Flexion	Posterior Offset	Modular	Force Options	No Grinding Necessary	Noise Dampening	Water Resistant