

# Leicht-Walker Walking Podium




## Features

Sizes available for the **smallest to large** patients

Lightweight **carbon construction**  
30% weight reduction

- + **Vertical alignment adjustment in four directions** without disassembly
- + **Parallel adjustment of step width** without disassembly
- + **Individual adjustment of step length** without disassembly
- + **Noise dampening movement stops**



## Standard

**Fully assembled** as one item for clear cost estimates and uncomplicated ordering

- + **Pyramid adapters** for precise adjustment
- + **Easy removal** of the body shell
- + **Height adjustable** to accommodate growth
- + **Secure 3-point attachment** of the body shell at the feet and pelvic area

**Upon request** an **extra-wide connecting plate** can be supplied for more abduction (aluminium)

- + **Universal connecting plate** is interchangeable with the Foot-Walker and Stehbrett

Unrivalled **20-month guarantee**


Engraved **serial number** for traceability

**Clearly illustrated** instructions for use



Article No.	Height Classification	Connecting Plate	Heel to Heel (middle)
Leicht-Walker-klein	60 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Leicht-Walker-mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Leicht-Walker-groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

### With anti-tip protection to prevent toppling backwards



Article No. Combi-product	Height Classification	Connecting Plate	Heel to Heel (middle)
Leicht-Walker- <b>KS</b> -klein	60 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Leicht-Walker- <b>KS</b> -mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Leicht-Walker- <b>KS</b> -groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

### With hand grips



Article No. Combi-product	Height Classification	Connecting Plate	Heel to Heel (middle)
Leicht-Walker- <b>HG</b> -klein	60 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Leicht-Walker- <b>HG</b> -mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Leicht-Walker- <b>HG</b> -groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

### With hand grips and anti-tip protection




Article No. Combi-product	Height Classification	Connecting Plate	Heel to Heel (middle)
Leicht-Walker- <b>HG-KS</b> -klein	60 cm to 100 cm	25 cm x 16 cm	23 cm to 8 cm
Leicht-Walker- <b>HG-KS</b> -mittel	80 cm to 130 cm	32 cm x 19 cm	30 cm to 12 cm
Leicht-Walker- <b>HG-KS</b> -groß	120 cm to 160 cm	39 cm x 22 cm	36 cm to 18 cm

### Body Shell with Connecting Plate

Custom body shell used in conjunction with the Stehbrett supports primary deformities and selectively counters any possible secondary deformities.

Only a precise footbed and proper support of the ankle joint enable correct positioning of the knee and hip joints.




Noise Dampening




Growable



Vertical Axis



Step Width



Step Length



3-Point Fixation



Removable



Reduced Weight



**HG** = Hand Grips



**KS** = Anti-tip Protection